CARDIOPULMONARY RESUSCITACION



- address and gently shake by the shoulders
 - perform a head tilt to open airway
 - 3. check the breathing no breathing? or abnormal breathing?
 - 4. call the ambulance
 - perform 30 chest compressions in the middle of the chest using both hands, depth of 5-6 cm, speed of 100-120 compressions/minute
 - deliver 2 rescue breaths
- 7. continue resuscitating in30:2 ratio until ambulance arrives











USE OF AED









AUTOMATED EXTERNAL DEFIBRILLATOR



- 2. start CPRimmediately
 - 3. if the AEDis present, take it (or send someone to bring it)
 - 4. open AEDand turn it on
 - ⁵ place electrodes as show on the pictograms
 - follow instructions given by AED





MASSIVE BLEEDING CONTROL





mind your own safety - use gloves or plastic bag as protection

find the source of bleeding

apply direct pressure on the wound with your finger or hand

4. if the first aid kit is available, apply pressure bandage

5. call the ambulance

monitor the wounded place until ambulance arrives







